

Delivering simplicity and affordability through innovation.

The SystemSURE Plus ATP luminometer is the next generation of the best selling SystemSURE II. Designed with state of the art electronics and improved functionality this palm-sized system is easy to use, extremely sensitive and affordable. Used by large and small food processors, hospitals, restaurants, supermarkets, and other industries where rapid detection of contamination is critical, the SystemSURE Plus allows companies to quickly determine the cleaning efficiency and hygienic status of surfaces and water, ensure product quality, and reduce costs.



How Does the SystemSURE Plus Work?

In conjunction with the Ultrasnap™ and Aquasnap™ testing devices, SystemSURE Plus measures adenosine triphosphate [ATP], the universal energy molecule found in all animal, plant, bacterial, yeast, and mold cells. Residues, particularly food or organic residue, contain large amounts of ATP. Microbial contamination also contains ATP, but in smaller amounts. After cleaning, all sources of ATP should be significantly reduced. When ATP is brought into contact with Hygiena's unique, liquid-stable luciferase/luciferin reagent in the Ultrasnap or Aquasnap testing device, light is emitted in direct proportion to the amount of ATP present. SystemSURE Plus measures the amount of light generated and provides information on the level of contamination in just seconds.


The higher the reading, the more contamination present. Assessing the cleanliness of a surface immediately after cleaning ensures contamination has been removed. Depending on the industry, this system will help product quality, extend product shelf life, prevent cross-contamination, enable immediate corrective action and/or avoid recalls. ATP hygiene monitoring provides accurate and traceable verification of the hygienic status of a surface or water, which is a key component of a good hygiene program. ATP testing is a universally recognized tool used by the world's largest and smallest companies for measuring the hygienic status of surfaces and water in order to ensure product consistency and safety.

SystemSURE Plus™

ATP Hygiene Monitoring System

Features:

- Small and lightweight
- Sensitive - detects down to 1 femtomole of ATP
- 250 programmable locations per sample plan
- 20 programmable sample plans
- Alphanumeric location ID
- Stores 2000 results
- Simple to use, menu-driven operation from keypad
- Results in 15 seconds
- Self-calibration at startup
- 2x AA battery powered - 3,000+ tests
- SURE-Trend analysis software (optional)



ATP Surface Test

Ultrasnap™



ATP Water Test

AQUASNAP™

New Technology Makes Ultrasnap and Aquasnap More Precise and Affordable

Ultrasnap and Aquasnap are user-friendly, self-contained ATP testing devices used with SystemSURE Plus. Both contain a unique liquid-stable reagent in place of the freeze-dried enzyme pellet used in other ATP testing devices. This new technology provides more accuracy, longer-lasting signal strength, and better reproducible results.

Features for Both:

- All-in-one sampling device
- 12-month shelf life
- Unique liquid-stable reagent
- Tolerant to temperature abuse
- Costs up to 50% less than other ATP tests




SURE-Trend ✓

Data Analysis Software


SURE-Trend is an easy-to-use result database management software. It allows companies to track results, quickly identify problem areas and generate reports for management or recordkeeping. Companies that implement SURE-Trend will reduce costs and increase productivity.

Features and Benefits


- ✓ PC/Windows® based software
- ✓ Intuitive user interface and operation
- ✓ Quick querying and trending of results
- ✓ Customize reports and charts
- ✓ Export data into Excel
- ✓ Identifies problem areas for immediate corrective action
- ✓ Saves time and lowers costs



Customize graphs to your needs



Quick statistics of location or sample plans



Easy set up of locations and pass/fail thresholds


SIMPLE
PRECISE
AFFORDABLE

4 easy steps | results in seconds



1

Swab



2

Snap & Squeeze



3

Insert



4

Read Results